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Routine Control of Steam Sterilizers in Hospitals by Means of Chemical Indicators

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Summary

Potential problems of pre-vacuum steam sterilization processes are discussed. The periodical control of the sterilization process by means of bio-indicators is necessary but not sufficient. For routine control, a Bowie-Dick-Test monitoring air removal and steam penetration and a batch control with chemical indicators are suggested as well as a simple documentation process of the retrieved information.

1 Introduction

The pre-vacuum sterilization process is the most commonly used procedure for the sterilization of multiple-use sterile goods. The advantage of the steam sterilization process is the use of inexpensive and non-toxic steam and, in comparison to all other sterilization types, extremely rapid sterilization cycles. As a result of this process, the specific investment costs, necessary man power and cost of steam is low. The only limiting factors are thermolabile products, products with corrosive surfaces and optical instruments. Consequently, more than 90% of all goods are sterilized in steam sterilization cycles.

In German speaking countries, the control of steam sterilizers with biological indicators is required between 1 and 6 months or 100 to 400 batches. The experience with vacuum steam sterilizers has shown that a periodical check with biological indicators is absolutely required, however, not sufficient to ensure safe routine operation of sterilizers in hospitals.

Modern fractionated vacuum sterilizers are able to sterilize goods at 134 °C between 4 and 7 minutes to 10^{-6} CFU per unit, if the bioburden of the sterilized goods is sufficiently low. Current wash and disinfection procedures guarantee such low bioburdens.

Short sterilization cycles between 4 and 7 minutes at 134 °C, however, are only sufficient, if the removal of air or non-condensable gases (e.g. CO₂) is ensured. In case of insufficient air removal, the sterilization cycle will not be safe even if all other parameters like steam pressure and temperature are in line with the given process parameters. Experience has shown that residual air or non-condensable gases in steam will be concentrated in linen packs and/or containers where high heat capacity is required. Due to the high consumption of steam to heat up the sterilization goods, large amounts are required inside these packages or containers. While the steam will condense on these goods, all other non-condensable gases like air in the steam will be transferred to the package and remain there without condensing. In these critical areas, the gas bubbles of non-condensable gases are isolating the goods to be sterilized and block further condensation of the steam. This will stop the rapid heat-up process, and in the gas bubbles, the heat-up phase becomes as long as in the well-known dry heat sterilizers.

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2 Potential Problems with Inadequate Air Removal

1. Failing vacuum pumps
Inadequate air removal of vacuum pumps increases the amount of residual air in the sterilization chamber. Even with vacuum pumps creating an evaporation pressure of 250–300 mbar, sufficient air removal may be reached, if 3–4 vacuum steps with sufficient steam pulses in between are programmed. These deficiencies may easily be recorded by existing pressure gauges. Therefore the risk is predictable.
2. Number of vacuum cycles
Depending on vacuum and chamber volume, a sufficient number of vacuum cycles is required. The number of cycles is usually determined by the sterilizer manufacturer with the initial validation of the system with biological indicators. There are a number of problems which may occur. For example, it happens very often that original standardized sterilization cycles are later modified without validating the total sterilization cycle again. Routine checks with biological indicators are usually carried out in standard cycles. Problems will not show up with modified programs.
3. Leakages of door seals, valves and other devices
A hole of 1 mm diameter is sufficient to negate the effect of a sterilization cycle. Door seals may leak due to various reasons: the surface is contaminated, the opposite metal surface is uneven, the seal broken, etc. Beginning problems cannot be easily detected. Only vacuum tests with an adequate vacuum holding time may detect small leaks. Such leak detect programs are seldom built in.
4. Pressure of air seal
The door is always sealed by air or steam pressure behind the seal. In principle, steam is to be preferred due to the fact that in case of leaks, steam causes no problems in relation to air penetration. However, using steam, seals have to withstand extremely high temperatures and will reduce the seal life. Door seals may not leak during a vacuum test, but this may be the case under steam pressure conditions. Also, leaks of valves behind seals may reduce the pressure and cause problems during the cycle not observed in the vacuum phase. This type of problem is not visible at the pressure gauges at all.
5. Non-condensable gases in steam
In hospitals, usually centralized process steam is used. This steam, however, is not adequately monitored for its content of inert gases. Reliable test procedures for steam quality are rarely available. Corresponding problems are often unknown. Inert gases in steam will heavily increase the total amount of non-condensable gases during the sterilization cycle. To minimize non-condensable gases, de-ionized water is prepared by ion exchange devices. If the ion exchanger loses its exchange capacity, hydrogen carbonate containing water will produce CO₂ and carbonate during steam production. This type of failure happens very often in hospitals. Sometimes the problems come only periodically. The daily Bowie-Dick-test carried out in the morning may be o. k. Failures later in the day remain undetected.

As yet there are no reliable measure devices available, indicating insufficient air removal. Monitoring temperature and pressure against time is important for sterilization process control, however, not sufficient to guarantee a safe sterilization cycle.

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3 Conclusions to be Drawn from Potential Malfunctions in the Vacuum Steam Sterilization Process

Referring to above mentioned problems, insufficient air removal, leaks or inert gases in the steam cause defective sterilization cycles. These problems cannot be detected with existing control devices (pressure and temperature gauges). Periodical checks with biological indicators only indicate problems if:

1. Problems are permanent and not periodical
2. If a permanent problem occurs after the prior check, this fault is not detected until the next periodical check is carried out.

To minimize the above mentioned risks, control of each individual sterilization cycle is required. Using biological indicators for that purpose is not practical because the sterilized items have to be stored in quarantine for several days. Therefore, other test methods have been established to show results immediately after the sterilization cycle:

4 Steam Penetration Test Detecting Non-Condensable Gases (Bowie-Dick-Test)

Bowie and Dick were the first ones to describe a textile pack of 6.6 kg with an enclosed indicator test sheet. This test is carried out after the start up phase in an empty sterilizer. It detects insufficient air removal and it is described in DIN 58946 part III. For practical use in the daily hospital environment, this test has some disadvantages:

1. Preparation of the test pack is labor-intensive and uneconomical
2. The reproducibility is influenced by:
 - Different packing density
 - Different age and porosity of textiles
 - Humidity of textiles
 - Residues of washing powder and softener ingredients of textiles
 - External wrapping
3. The test results are only valid, if the sterilization cycle at 134 °C of 3.5 min is not exceeded. Therefore sterilizers require a special BD-test program. Hospitals which do not have this BD-test program cannot use the original BD-test. Running the BD-test in a longer cycle does not prove anything.

To circumvent the production of individual Bowie-Dick-test packs, various Bowie-Dick-Simulation tests (BDS-tests) have been developed, giving acceptable results in most cases. Unfortunately, today necessary test methods for these BDS-tests are not available. Single-use BDS-tests are fairly well reproducible. The cost per test (DM 5,- to DM 10,-) is relatively high. Moreover, paper or plastic waste is produced. Multiple-use BDS-tests (usually two paper or non-woven packs with a Bowie-Dick-test sheet in between) are more economical. The disadvantage of these packs is that the porosity of the packs is modified during each sterilization cycle. Therefore, the reproducibility is slightly modified by each cycle. These packs are exchanged after 15–25 cycles.

GKE developed a test device based on a metal carrier and a plastic helix which can be used a number of times without any change of the test device itself. The long-term reproducibility of this inexpensive test device is guaranteed. Waste is not produced. The documentation of a fairly small indicator test strip is easy.

The above mentioned Bowie-Dick-test does not represent a sterility test for the sterilizer. This test only indicates that all

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physical requirements in the vacuum steam sterilizer are obtained to achieve a sterilizing cycle, killing all germs in a short period of time. The daily Bowie-Dick-test increases the probability of acceptable sterilization cycles. However, a 100% probability cannot be given, as the above described malfunctions may occur after the test.

5 In-line Batch or Cycle Control

The individual control of each load considerably increases the safety of a sterilization process. Only the simultaneous break down of sterilizer and cycle control system would result in an undetected faulty sterilization cycle. To use biological indicators for cycle control is impractical as described above. In 1986, GKE developed a cycle control system with chemical indicators being able to check faulty sterilization processes after validation with biological indicators. The advantage of chemical indicators getting the results immediately after the sterilization cycle is obvious.

Very often self-adhesive indicators tapes, indicators printed on pouches, or chemical indicators attached to containers are confused with cycle control systems. This type of indicators are process indicators, exclusively giving the information that those packs/containers have passed the sterilization process without any indication of the process quality or sterility. Many operating theatre nurses believe that a changed chemical indicator colour on the surface of a package or container is equal to a sterile contents. In reality, it is not a sterility check. The only information given is that the product was subjected to a processing cycle. Even integrative indicators do not give sterility assurance, if they are fixed on the surface of packages, containers, or if the indicator itself is placed between packages. Chemical cycle control systems are only able to check the quality of a sterilization process successfully, if the following conditions are met:

1. The chemical indicator should change color only if the necessary sterilization conditions are achieved, i.e. >134 °C, 3 min, 121 °C 15 min.
2. The conditions of point 1 need to be achieved in areas where worst case sterilization conditions occur. If air is present in these areas, sterility conditions are not achieved or the exposure time is not sufficient to transfer the necessary heat capacity in densely packed containers or packs. Therefore the placing of a chemical indicator requires:
 - a densely packed package or container or
 - using a test device, containing the chemical indicator and simulating worst case sterilization conditions outside of a pack.

For that purpose, temperature-time-integrating indicators are used. They are placed in the worst case location of a textile pack or container. The disadvantage of this method is that the result of the sterility test will be available only after the package or container is opened for use.

To overcome this disadvantage, simulation test devices with incorporated time-temperature-integrating chemical indicators have been developed to simulate worst case conditions. The combination of indicator and test device is called batch or cycle control test device. These systems are validated against biological indicators. The in-line batch or cycle control systems can be used as worst case indicators to check the whole contents of a batch.

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6 Control of all Used Sterile Medical Devices

Sterile medical devices are produced under GMP (Good Manufacturing Practice) or CEN EN 29 XXX and EN 45 XXX rules in medical industry. These rules give a detailed description of production and validation control. Validated procedures require continuing quality control. Companies working under these conditions are able offer the user a quality certificate with the guarantee of sterility. Hospitals having quality certificates can use medical devices without documenting all individual products used.

Most hospitals do not sterilize medical products under GMP rules. Therefore, general quality certificates cannot be passed to the end-user. However, similar requirements as defined under GMP must be implemented:

1. Definition of information flow
2. Definition of materials flow
3. Definition of responsibilities
4. Definition of quality control
5. Definition of documentation
6. Definition of the production process
7. Validation of processes
8. Routine control

It is obvious that under hospital conditions, GMP-rules cannot be easily implemented. However, the goal to produce exclusively sterile products is the same as in the medical industry.

7 Documentation of the Results after Bowie-Dick and Cycle Control Tests

For documentation of the Bowie-Dick- and cycle control-tests, load record sheets or cards are offered to collect all information of one day and sterilizer for the Bowie-Dick and all cycle control tests. The indicator strips are self-adhesive and placed on the record card. The person in charge checks the results and writes them on this documentation sheet. All packs or containers of a cycle are labeled with the production date, sterilizer and batch number. This label is only placed on the pack or container if the cycle control system indicates sterility conditions. This label signals sterility to the end-user. The use of process indicators on each individual package is not required.

The quick production and transfer of these labels is achieved by labeling print devices similar to those used in supermarkets.

If a patient orientated documentation of all individual sterile goods is desired, double self-adhesive labels can be used. Before opening the package in the operation theatre, the label can be transferred to the patient documentation file. In this case, the sterilized items can be traced back to the respective load record sheets. ■

Dr. Ulrich Kaiser
GKE mbH
Feldstrasse 14
D-61479 Glashütten
Germany