Containerplombe | Container-Seal

Art.-No.	Content [pcs]	Product Code	Product Description	Application
211-362	1.008	C-S-L-1-CS	self-adhesive container seal with process indicator for steam sterilization processes 29 x 70 mm, pin length 75 x 3,5 mm	to seal metal containers with small holes for the seals

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The container seal is made from self-adhesive carton with a dimension of 7 x 2,9 cm and a pin length of 7,5 cm. The pin is passed through the seal mechanism of the container and is fixed between two self-adhesive layers. The adhesion is steam resistant, so that the container cannot be opened without destroying the seal.

The seal contains a steam process indicator that a destroyed seal cannot be exchanged without notice. Also, logistic information is available, which containers have been processed already. The seal label has enough space to write all necessary information to fulfil all QC requirements.

Alternatively, the GKE hand labelling device, art.-no. 240-820 or 240-830, may be used in combination with the labels 240-861 to 240-864 to replace the hand-written information. After usage on the container in the operation room the double self-adhesive label may be taken off and adhered to the patient documentation sheet.

The described procedure may replace conventional container marking. The paper-based container seals are easily recycled together with conventional paper. No specific waste collection is required.

760-021 ED V02