

Practise, Laboratory, Hospital

.....

.....

Documentation Sheet

for incubation of GKE Steri-Record® Instant-Mini-Bio-Plus
self-contained biological indicators (SCBI)
for steam sterilization processes



Batch no. of SCBI		Result vitality test (***)	STERILE		NON-STERILE	
Expiry Date of SCBI		Incubation temperature	55 - 60°C			
SCBI no. for vitality test*		Maximum recommended incubation time	24 hours			

No. (*)	Date of sterilization	Sterilizer No.	Batch No. of sterilization process tested	Type 5 indicator (**) Read-out result		Incubation period				Result after incubation (***)		Responsible Person	
				Sterile	Non-sterile	Start		End		Sterile	Non-sterile	Initials	Signature
						Date	Time	Date	Time				

766-002 EN V03 03/2023

(*) If several Instant-Mini-Bio-Plus SCBIs and/or vitality test are incubated as well, all ampoules shall be marked with a number before use, to assign the results correctly. Insert this number in this column.

Explanation: Please tick-off:

- (**) STERILE = All 3 indicator bars have changed from yellow to black.
NIGHT STERIL = Indicator bars have not or not fully changed from yellow to black.
- (***) STERILE = No growth, no colour change of the growth medium
NON-STERILE = Growth with colour change of the growth medium