

# Documentation Sheet



## GKE Clean-Record® CPI for bedpan washers

Hospital/Clinic:	City:
Department:	Date:
Manufacturer:	Type:
Year of manufacture:	Machine No.:
Indicator Lot:	Checked by:

The below documented test result is

☐ a **benchmark result** for the above-mentioned bedpan-washer. Later on, tests are comparable in the identically constructed machines.

☐ daily ☐ weekly **routine testing**

**Evaluation** ☐ Test result = reference → OK

☐ Test result ≠ reference → not OK

Program:

Lot-No.:

Bedpan					
Bottom, inside, innen	On the edge, right	On the edge, left	Under the edge, right	Under the edge, left	Bottom, outside
<i>Adhere indicator here</i>					
<b>Handle</b>	<i>Additional test position</i>	<i>Additional test position</i>	<i>Additional test position</i>	<i>Additional test position</i>	<i>Additional test position</i>

Urine bottle					
Lateral, outside	Bottom, outside	Bottle neck, inside	Lateral, inside	Bottom inside	<i>Additional test position</i>
<i>Adhere indicator here</i>					