

**Documentation Sheet for the results of
GKE Steri-Record® Testset (200-013 ED)
to check the penetration characteristics of sterilization processes**

Hospital/ Clinic: _____ Sterilizer: _____

Address: _____ Type: _____

_____ Date of manufacture: _____

Program: _____ Date: _____

Air removal:

Sterilization process: ☐ steam ☐ formaldehyde ☐ hydrogen peroxide

Sterilization: Temp. °C Holding time min

PCD-No.	Art.-No.	Tube Dimensions		HPR* [cm²]	Colour Change of chemical indicator
		Length [cm]	Diameter [cm]		
1	200-150	150	0,2 (= 2 mm)	30	Please adhere chemical indicator here
2	200-100	100		20	Please adhere chemical indicator here
3	200-075	75		15	Please adhere chemical indicator here
4	200-050	50		10	Please adhere chemical indicator here
5	200-025	25		5	Please adhere chemical indicator here

Program: _____ Date: _____

Air removal:

Sterilization process: ☐ steam ☐ formaldehyde ☐ hydrogen peroxide

Sterilization: Temp. °C Holding time min

PCD-No.	Art.-No.	Tube Dimensions		HPR* [cm²]	Colour Change of chemical indicator
		Length [cm]	Diameter [cm]		
1	200-150	150	0,2 (= 2 mm)	30	Please adhere chemical indicator here
2	200-100	100		20	Please adhere chemical indicator here
3	200-075	75		15	Please adhere chemical indicator here
4	200-050	50		10	Please adhere chemical indicator here
5	200-025	25		5	Please adhere chemical indicator here

*HPR [cm²] = Length [cm] x Diameter [cm]